


# Helpful Tips for Laboratories

**Thank you for choosing Anatotemp SC®, the first and only anatomic dental implant healing abutment and scanbody. Anatotemp SC® is an off-the-shelf solution for creating ideal emergence profile and streamlines the digital dental implant workflow.**

The Anatotemp SC is used like any other digital scan body/post. The clinician scans the case and sends scan along with the Anatotemp SC identification number as well as implant identification and platform size or nomenclature.

As a reminder, Anatotemp SC comes in six distinct shapes and multiple internal anti-rotational connections as described below:



Laboratories must download the most recent CAD and Model Creator libraries for their software. (i.e. Exocad, 3Shape)

Laboratories must choose the appropriate implant analog based on implant brand/connection and corresponding analog.

#### EXAMPLES:

- Internal Hex Connection → Biohorizons implant analog ( PYIA, PGIA, PBIA)
- Conical Hex Connection → Nobel Active implant analog (38189, 38190, 38191)
- Internal 4Side Connection → Straumann Crossfit implant analog (0254101)
- Ditron Hex Connection → Ditron implant analog (IDA - 2030)

Often, one must validate the Analog - Model fit by adjusting the offset (Exocad) or friction bars (3Shape) when utilizing the model creating software.

