



INTERNAL OCTAGON IMPLANT CHART

4.8mmD INTERNAL OCTAGON

IMPLANT COMPANY	IMPLANTS
AS Technology (Titanium Fix)	c-fix Cone Morse Regular (4.8)
B&W s.r.l.	Tapered Morse MEI (4.8)
Biomain AB	I-Fix Regular (4.8)
Biomet 3i	TG OSSEOTITE (4.8)
Blue Sky Bio	One Stage Narrow Platform (4.1), One Stage Regular Platform (4.8)
Clinical House Europe GmbH	PerioType TG (4.8)
Cowellmedi Corporation, Ltd	Atlas Internal Regular, Atlas Internal Taper (All 4.8)
CSM	Apolonia Internal Straight Type, Apolonia Internal Tapered Type (All 4.8)
Dental Ratio Systems (DRS) GmbH	Oktagon (4.8)
Dentis	Internal
Dr. Ihde Dental GmbH	Allfit OSS, Allfit Osseopore, Allfit SSO, Allfit STC, Allfit STI, Allfit STI RT, Allfit STO, Allfit TPG (All 4.8)
Euroteknika	Aesthetica (4.2), Aesthetica (4.8), Aesthetica+ (NP-4.2), Aesthetica+ (RP-4.8)
Hi-Tec Implants	Tapered Tite Fit, Tite Fit (All 4.8)
Implant Direct LLC	SwishPlant (4.8)
Intra-Lock System International	Standard Platform Conic Body (4.3), Standard Platform Straight Body (4.3), Wide Platform Conic Body (5.5), Wide Platform Straight Body (5.5)
Keystone Dental (Lifecore)	Stage-1 RDS
MegaGen	ExFeel Internal Straight (4.8), ExFeel Internal Taper (4.8)
Microtech srl	Microfit (4.8)
MIS Implant Technologies	MISTRAL (4.8)
Neobiotech Co., Ltd	Neoplant Internal (4.8)
Osstem Implant Co., Ltd (Hiossen)	SS I Regular, SS II Regular, SS III Regular (All 4.8)
Renew Biocare AG	Re-Fit TG Octa (4.8), Re-Fit TG ST (4.1)
Southern Implants	ITC, ITSA, ITT, Max IT 7, OCT ITSTA, OCT-XT ITST 12d Co-Axis (All 4.8)
Sterngold Dental, LLC	Stern IC (4.8)
Straumann AG	Hollow Cylinder 15 Degree Angled Esthetic Plus 3.5 (4.8), Hollow Cylinder Esthetic Plus 3.5 (4.8), Standard (RN), Standard Plus (RN), Taper Effect (RN)
Sweden & Martina SPA	Stark-D (4.8)
Sybron Implant Solutions	Endopore One Stage, SybronPRO TL, SybronPRO XRT (All 4.8)
Trinon Titanium	QK
UFIT Implant Inc.	ST1, VT1 (All 4.8)
Zimmer Dental	SwissPlus Straight, SwissPlus Tapered (All 4.8)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.



INTERNAL OCTAGON IMPLANT CHART

6.5mmD INTERNAL OCTAGON

IMPLANT COMPANY	IMPLANTS
AS Technology (Titanium Fix)	
B&W s.r.l.	
Biomain AB	
Biomet 3i	TG OSSEOTITE (6.5)
Blue Sky Bio	One Stage Wide Platform (6.5)
Clinical House Europe GmbH	
Cowellmedi Corporation, Ltd	
CSM	
Dental Ratio Systems (DRS) GmbH	
Dentis	
Dr. Ihde Dental GmbH	Allfit STW (6.5)
Euroteknika	Aesthetica (6.5), Aesthetica+ (WP-6.5)
Hi-Tec Implants	
Implant Direct LLC	SwishPlant (6.5)
Intra-Lock System International	
Keystone Dental (Lifecore)	Stage-1 WDS
MegaGen	
Microtech srl	
MIS Implant Technologies	
Neobiotech Co., Ltd	
Osstem Implant Co., Ltd (Hiossen)	
Renew Biocare AG	
Southern Implants	ITC 6 (6.5), ITT 6 (6.5), Max IT 8 (6.5), Max IT 9 (6.5,7.5), OCT ITST6 (6.5)
Sterngold Dental, LLC	
Straumann AG	Standard (WN), Standard Plus (WN), Taper Effect (WN),
Sweden & Martina SPA	
Sybron Implant Solutions	Sybron PRO TL (6.5), SybronPRO XRT (6.5)
Trinon Titanium	
UFIT Implant Inc.	ST1 (6.5), VT1 (6.5)
Zimmer Dental	

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.