


INTERNAL TRI-LOBE IMPLANT CHART

3.5mmD INTERNAL TRI-LOBE

IMPLANT COMPANY	IMPLANTS
ACE Surgical	infinity TRI-CAM (3.5)
Biomain AB	Combine Narrow (3.5)
Blue Sky Bio	Two Stage Trilobe (3.5)
Hi-Tec Implants	Vision VRI (3.5)
Implant Direct LLC	ReActive, RePlant, RePlus (All 3.5)
Nobel Biocare AB	Nobel Replace - Straight NP, Nobel Replace - Taper NP, Nobel Speedy Replace (All 3.5)
Renew Biocare AG.	Re-Root II SD (3.5)
Southern Implants	Tri-Nex Co-Axis IA43-12d (3.5), Tri-Nex IA-L-35, Tri-Nex IA-LHS-35
Tecom Implantology (TITANMED S.r.l.)	REP STR GR TU (3.5), REP Tapered GR TU (3.5)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.





INTERNAL TRI-LOBE IMPLANT CHART

4.3mmD INTERNAL TRI-LOBE

IMPLANT COMPANY	IMPLANTS
ACE Surgical	infinity TRI-CAM (4.3)
Biomain AB	Combine Regular (4.3)
Blue Sky Bio	Two Stage Trilobe (4.3)
Hi-Tec Implants	Vision VRI 4.3
Implant Direct LLC	ReActive, RePlant, RePlus (All 4.3)
Nobel Biocare AB	Nobel Replace - Straight RP (4.3), Nobel Replace - Taper RP (4.3), Nobel Speedy Replace (4.3)
Renew Biocare AG.	Re-Root II MD (4.3)
Southern Implants	Tri-Nex Co-Axis IA50-12d (4.3), Tri-Nex IA-L-43, Tri-Nex IA-LHS-43
Tecom Implantology (TITANMED S.r.l.)	REP STR GR TU (4.3), REP Tapered GR TU (4.3)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.





INTERNAL TRI-LOBE IMPLANT CHART

5.0mmD INTERNAL TRI-LOBE

IMPLANT COMPANY	IMPLANTS
ACE Surgical	infinity TRI-CAM (5.0)
Biomain AB	Combine Wide (5.0)
Blue Sky Bio	Two Stage Trilobe (5.0)
Hi-Tec Implants	Vision VRI 5.0
Implant Direct LLC	ReActive, RePlant, RePlus (All 5.0)
Nobel Biocare AB	Nobel Replace - Straight WP, Nobel Replace - Taper WP, Nobel Speedy Replace (All 5.0)
Renew Biocare AG.	Re-Root II WD (5.0)
Southern Implants	Tri-Max 7 (5.0), Tri-Nex IA-L-50
Tecom Implantology (TITANMED S.r.l.)	REP STR GR TU (5.0), REP Tapered GR TU (5.0)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.



INTERNAL TRI-LOBE IMPLANT CHART

6.0mmD INTERNAL TRI-LOBE

IMPLANT COMPANY	IMPLANTS
ACE Surgical	
Biomain AB	
Blue Sky Bio	
Hi-Tec Implants	
Implant Direct LLC	RePlant (6.0)
Nobel Biocare AB	Nobel Replace - Taper 6.0 (6.0), Nobel Speedy Replace (6.0)
Renew Biocare AG.	Re-Root II WD (6.0)
Southern Implants	Tri-Max 8 (6.0), Tri-Max 9 (6.0), Tri-Nex IA-L-60
Tecom Implantology (TITANMED S.r.l.)	

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.