

# Anatotemp Dimensions

Shape	Height*	Width (Buccal-Lingual)	Width (Mesial-Distal)
Maxillary Central	5 mm	7.5 mm	7 mm
Maxillary Lateral - Mandibular Incisor	5 mm	6.5 mm	5 mm
Maxillary Canine	5 mm	8.25 mm	6.5 mm
Maxillary/Mandibular Premolar & Mandibular Canine	5 mm	8.5 mm	5.75 mm
Maxillary Molar	5 mm	9 mm	9 mm
Mandibular Molar	5 mm	8.5 mm	9.5 mm

\*The **Height** dimension is the height of the Anatotemp when fully seated on the implant.