

# AnatotempSC<sup>®</sup>

Anatomic Dental Implant Healing Abutment and Scan Body

Three dots on the buccal!



The Anatotemp anatomic healing abutment has become a staple in providing an off-the-shelf solution for creating ideal gingival emergence profile. We now introduce AnatotempSC, the second generation Anatotemp that is both an anatomic healing abutment and a digital impression body. The digital dental implant workflow has become even easier saving four steps and one appointment by utilizing the AnatotempSC. With the AnatotempSC there is no need to remove the healing abutment and place an impression post or scan body. AnatotempSC does it all! AnatotempSC works with many major dental implant systems.

## Six Anatomic Shapes



## Many Anti-Rotational Connections




The off-the-shelf solution for creating ideal emergence profile!

The AnatotempSC can be ordered through your Ditrion Dental implant representative.



# AnatotempSC® - Shape Location Guide

The Anatotemp Shape Location Guide provides guidance for identifying which Anatotemp shape to utilize for each tooth replaced.



## AnatotempSC - Product Selection

The AnatotempSC - Ditron Internal Connection comes in six anatomical shapes. The matrix below organizes the Ditron AnatotempSC product numbers with the appropriate Ditron ULT & MPI implants.

### DITRON INTERNAL HEXAGON CONNECTION



	ANATOTEMP SC SHAPE					
	Mandibular Molar	Maxillary Molar	Maxillary/Mandibular Premolar Mandibular Canine	Maxillary Canine	Maxillary Lateral/Mandibular Incisor	Maxillary Central
ULT & MPI						
Diameter						
3.3, 3.5			DDI - 04SC	DDI - 03SC	DDI - 02SC	DDI - 01SC
3.5, 3.7, 4.2			DDI - 04SC	DDI - 03SC		DDI - 01SC
5.0, 6.0	DDI - 06SC	DDI - 05SC				

Compatible with: Ditron Dental USA ULT & MPI Implants Systems