



# AnatotempSC® TiBase - Shape Location Guide

The AnatotempSC Shape Location Guide provides guidance for identifying which Anatotemp shape to utilize for each tooth replaced.



## Tooth Numbers



# AnatotempSC® TiBase - Dimensions

Shape	Height	Width (Buccal-Lingual)	Width (Mesial-Distal)
Maxillary Central	5 mm	7.5 mm	7 mm
Maxillary Lateral - Mandibular Incisor	5 mm	6.5 mm	5 mm
Maxillary Canine	5 mm	8.25 mm	6.5 mm
Maxillary/Mandibular Premolar & Mandibular Canine	5 mm	8.5 mm	5.75 mm
Maxillary Molar	5 mm	9 mm	9 mm
Mandibular Molar	5 mm	8.5 mm	9.5 mm

Anatotemp shapes and sizes were designed utilizing an exhaustive cone beam CT study, enabling the formulation of average sizes that are ideal for an off-the-shelf component.