


INTERNAL CONICAL HEXAGON IMPLANT CHART


3.0mmD INTERNAL CONICAL HEXAGON

| IMPLANT COMPANY | IMPLANTS |
|--------------------|---|
| Nobel Biocare AB | Nobel Active-Internal Connection NP (3.5) |
| Blue Sky Bio | Max |
| Hi-Tec Implants | Logic LGI 3.5 |
| Implant Direct LLC | InterActive 3.0mmD Platform (3.2, 3.7), SwishActive 3.0mmD Platform (3.3) |

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.



Anatotemp[®]



INTERNAL CONICAL HEXAGON IMPLANT CHART

3.4mmD INTERNAL CONICAL HEXAGON

| IMPLANT COMPANY | IMPLANTS |
|--------------------|--|
| Nobel Biocare AB | Nobel Active - Internal Connection RP (4.3, 5.0) |
| Hi-Tec Implants | Logic LGI (4.3, 5.0) |
| Implant Direct LLC | InterActive Implant 3.4mmD Platform (4.3, 5.0), SwishActive 3.4mmD Platform (4.1, 4.8) |

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.