



INTERNAL HEXAGON IMPLANT CHART

3.5mmD INTERNAL HEXAGON

IMPLANT COMPANY	IMPLANTS
A.B. Dental Devices	I10, I15, I2 Narrow Platform, I2b Narrow Platform, I5 Narrow Platform (All 3.75)
ACE Surgical	infinity Internal Hex (3.7 & 4.1) (3.5mmD Platform)
ADIN Dental Implants System	Swell, Touareg, Touareg-s
Alfa Gate	Bioactive, DFCI (3.5)
Alpha Bio Tec.	ATID, DFI, ICE (3.5), SFB, SPI
Alpha Dent Implants	Alpha Dent Active, Active Bio, Active Conus, Active Plus, Classic, Classic Conus
American Dental Implant Corporation - ADI	Internal Hex Skinny 2.4, Internal Hex Small Diameter, Internal Hex Standard Diameter, Internal Hex Tapered Standard Diameter (All 3.5)
ARDS Implants	Classic, Dental Implant S, Premium, Smart (All 3.75)
BioHorizon	Internal, Single Stage, Tapered Internal, Tapered Internal Plus (All 3.5)
Blue Sky Bio	Internal Hex (3.5)
CORTEX Dental Implants Industries, Ltd.	Classic, Classix (3.75), Dynamic, Dynamix (3.75)
Den Tack	Quad Internal Hex
DENTIN Implants Technologies LTD	Classic, Prestige, Rapid
Dr. Ihde Dental GmbH	Allfit ATI(3.5)
Easy System Implant	Previum-S (4.1)
Glidewell Laboratories	Inclusive Tapered (3.5)
Hi-Tec Implants	Logic Plus, Smooth & Slim Fit, Tapered Self Thread, X6 (All 3.5)
Implant Direct LLC	Legacy 1, Legacy 2, Legacy 3, Legacy 4 (All 3.5)
Implasa-Höchst GmbH	DoubleGotic, Gotic, Nordic (All 3.5)
IQ Implants	Master, Smart (All 3.5)
Maxillent	iRaise, iSure (All 3.5)
Microtech srl	Minivent (3.5)
MIS Implant Technologies	BioCom Internal Hexagon - Standard Platform, Lance Internal Hexagon Standard Platform, Lance Plus - Standard Platform, M4 Standard Platform, Seven - Standard Platform (All 3.75)
Nova-Implants Ltd	Pacific, Spirit (All 3.5)
Park Dental Research Corp.	Star Vent Internal Hex Screw
Sialo Technology Israel Ltd.	DIVA (Dynamic Implant Valve Approach)
Sigma Implants Ltd	Sigma Active, Sigma Classic
TAV Dental	Silhouette, Sirius (All 3.5)
Tecom Implantology (TITANMED S.r.l)	Screw Straight, Screw Tapered, Starfly 0.2 Two Phases (All 3.5)
Tidal Spiral Implant Systems (TSI)	Tidal Screw Implant System TSS, Tidal Spiral Implant System TSI (All 3.5)
Top Implants LTD.	Dolphin, Orca, T-Shark (All 3.5)
TRI Dental Implants GmbH	TRI-Vent (3.5)
Unidental Direct, LLC	IN-SIX
Zimmer Dental	[Paragon] - Bio-Vent SD (3.5), Core-Vent SD (3.5), Micro-Vent SD (3.25), Micro-Vent Tapered SD (3.7) - [Paragon] AND Screw Vent, Tapered Screw Vent, Trabecular Metal (All 3.5)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.





INTERNAL HEXAGON IMPLANT CHART

4.5mmD INTERNAL HEXAGON

IMPLANT COMPANY	IMPLANTS
A.B. Dental Devices	I1 Wide Platform, I2 Wide Platform, I2b Wide Platform, I5 Wide Platform (All 4.5)
ACE Surgical	infinity Internal Hex (4.7 & 5.1) (4.5 mmD Platform)
ADIN Dental Implants System	
Alfa Gate	
Alpha Bio Tec.	
Alpha Dent Implants	
American Dental Implant Corporation - ADI	Internal Hex Tapered Wide Diameter, Internal Hex Wide Diameter (All 4.5)
ARDS Implants	Dental Implant S (4.5)
BioHorizon	Internal, Single Stage, Tapered Internal, Tapered Internal Plus (All 4.5)
Blue Sky Bio	Internal Hex, Molar (All 4.5)
CORTEX Dental Implants Industries, Ltd.	
Den Tack	
DENTIN Implants Technologies LTD	
Dr. Ihde Dental GmbH	
Easy System Implant	Previum-S (5.0)
Glidewell Laboratories	Inclusive Tapered (4.5)
Hi-Tec Implants	Tapered Self Thread (4.5)
Implant Direct LLC	Legacy 1, Legacy 2, Legacy 3, Legacy 4 (All 4.5)
Implasa-Höchst GmbH	
IQ Implants	
Maxillent	
Microtech srl	Minivent (4.7)
MIS Implant Technologies	BioCom Internal Hexagon Wide Platform, Lance Internal Hexagon Wide Platform, Lance Plus Wide Platform, M4 Wide Platform, One Stage Internal Hexagon, Seven Wide Platform (All 4.5)
Nova-Implants Ltd	Pacific (4.5)
Park Dental Research Corp.	
Sialo Technology Israel Ltd.	
Sigma Implants Ltd	
TAV Dental	
Tecom Implantology (TITANMED S.r.l)	Screw Straight (4.5)
Tidal Spiral Implant Systems (TSI)	Tidal Screw Implant System TSS, Tidal Spiral Implant System TSI (All 4.5)
Top Implants LTD.	Dolphin, Orca, T-Shark (All 4.5)
TRI Dental Implants GmbH	
Unidental Direct, LLC	
Zimmer Dental	[Paragon] - Bio-Vent MD (4.5), Core-Vent MD (4.5), Micro-Vent MD (4.25,4.7), Micro-Vent Tapered MD (4.7) AND Screw-Vent, Tapered Screw-Vent (All 4.5)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.




INTERNAL HEXAGON IMPLANT CHART

5.7mmD INTERNAL HEXAGON

IMPLANT COMPANY	IMPLANTS
A.B. Dental Devices	
ACE Surgical	
ADIN Dental Implants System	
Alfa Gate	
Alpha Bio Tec.	
Alpha Dent Implants	
American Dental Implant Corporation - ADI	Internal Hex Tapered XL (5.7)
ARDS Implants	
BioHorizon	Internal, Single Stage, Tapered Internal (All 5.7)
Blue Sky Bio	
CORTEX Dental Implants Industries, Ltd.	
Den Tack	
DENTIN Implants Technologies LTD	
Dr. Ihde Dental GmbH	
Easy System Implant	
Glidewell Laboratories	
Hi-Tec Implants	
Implant Direct LLC	Legacy 1, Legacy 2, Legacy 3, Legacy 4 (All 5.7)
Implasa-Höchst GmbH	
IQ Implants	
Maxillent	
Microtech srl	Minivent (5.7)
MIS Implant Technologies	
Nova-Implants Ltd	
Park Dental Research Corp.	
Sialo Technology Israel Ltd.	
Sigma Implants Ltd	
TAV Dental	
Tecom Implantology (TITANMED S.r.l)	
Tidal Spiral Implant Systems (TSI)	
Top Implants LTD.	
TRI Dental Implants GmbH	
Unidental Direct, LLC	
Zimmer Dental	[Paragon] - Micro-Vent Tapered WD (5.7) AND AdVent, Tapered Screw-Vent (5.7)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.




INTERNAL CONICAL HEXAGON IMPLANT CHART


3.0mmD INTERNAL CONICAL HEXAGON

IMPLANT COMPANY	IMPLANTS
Nobel Biocare AB	Nobel Active-Internal Connection NP (3.5)
Blue Sky Bio	Max
Hi-Tec Implants	Logic LGI 3.5
Implant Direct LLC	InterActive 3.0mmD Platform (3.2, 3.7), SwishActive 3.0mmD Platform (3.3)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.



Anatotemp[®]





INTERNAL CONICAL HEXAGON IMPLANT CHART

3.4mmD INTERNAL CONICAL HEXAGON

IMPLANT COMPANY	IMPLANTS
Nobel Biocare AB	Nobel Active - Internal Connection RP (4.3, 5.0)
Hi-Tec Implants	Logic LGI (4.3, 5.0)
Implant Direct LLC	InterActive Implant 3.4mmD Platform (4.3, 5.0), SwishActive 3.4mmD Platform (4.1, 4.8)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.





INTERNAL TRI-LOBE IMPLANT CHART

3.5mmD INTERNAL TRI-LOBE

IMPLANT COMPANY	IMPLANTS
ACE Surgical	infinity TRI-CAM (3.5)
Biomain AB	Combine Narrow (3.5)
Blue Sky Bio	Two Stage Trilobe (3.5)
Hi-Tec Implants	Vision VRI (3.5)
Implant Direct LLC	ReActive, RePlant, RePlus (All 3.5)
Nobel Biocare AB	Nobel Replace - Straight NP, Nobel Replace - Taper NP, Nobel Speedy Replace (All 3.5)
Renew Biocare AG.	Re-Root II SD (3.5)
Southern Implants	Tri-Nex Co-Axis IA43-12d (3.5), Tri-Nex IA-L-35, Tri-Nex IA-LHS-35
Tecom Implantology (TITANMED S.r.l.)	REP STR GR TU (3.5), REP Tapered GR TU (3.5)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.




INTERNAL TRI-LOBE IMPLANT CHART

4.3mmD INTERNAL TRI-LOBE

IMPLANT COMPANY	IMPLANTS
ACE Surgical	infinity TRI-CAM (4.3)
Biomain AB	Combine Regular (4.3)
Blue Sky Bio	Two Stage Trilobe (4.3)
Hi-Tec Implants	Vision VRI 4.3
Implant Direct LLC	ReActive, RePlant, RePlus (All 4.3)
Nobel Biocare AB	Nobel Replace - Straight RP (4.3), Nobel Replace - Taper RP (4.3), Nobel Speedy Replace (4.3)
Renew Biocare AG.	Re-Root II MD (4.3)
Southern Implants	Tri-Nex Co-Axis IA50-12d (4.3), Tri-Nex IA-L-43, Tri-Nex IA-LHS-43
Tecom Implantology (TITANMED S.r.l.)	REP STR GR TU (4.3), REP Tapered GR TU (4.3)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.



INTERNAL TRI-LOBE IMPLANT CHART


5.0mmD INTERNAL TRI-LOBE

IMPLANT COMPANY	IMPLANTS
ACE Surgical	infinity TRI-CAM (5.0)
Biomain AB	Combine Wide (5.0)
Blue Sky Bio	Two Stage Trilobe (5.0)
Hi-Tec Implants	Vision VRI 5.0
Implant Direct LLC	ReActive, RePlant, RePlus (All 5.0)
Nobel Biocare AB	Nobel Replace - Straight WP, Nobel Replace - Taper WP, Nobel Speedy Replace (All 5.0)
Renew Biocare AG.	Re-Root II WD (5.0)
Southern Implants	Tri-Max 7 (5.0), Tri-Nex IA-L-50
Tecom Implantology (TITANMED S.r.l.)	REP STR GR TU (5.0), REP Tapered GR TU (5.0)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.



Anatotemp[®]





INTERNAL TRI-LOBE IMPLANT CHART

6.0mmD INTERNAL TRI-LOBE

IMPLANT COMPANY	IMPLANTS
ACE Surgical	
Biomain AB	
Blue Sky Bio	
Hi-Tec Implants	
Implant Direct LLC	RePlant (6.0)
Nobel Biocare AB	Nobel Replace - Taper 6.0 (6.0), Nobel Speedy Replace (6.0)
Renew Biocare AG.	Re-Root II WD (6.0)
Southern Implants	Tri-Max 8 (6.0), Tri-Max 9 (6.0), Tri-Nex IA-L-60
Tecom Implantology (TITANMED S.r.l.)	

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.





INTERNAL OCTAGON IMPLANT CHART

4.8mmD INTERNAL OCTAGON

IMPLANT COMPANY	IMPLANTS
AS Technology (Titanium Fix)	c-fix Cone Morse Regular (4.8)
B&W s.r.l.	Tapered Morse MEI (4.8)
Biomain AB	I-Fix Regular (4.8)
Biomet 3i	TG OSSEOTITE (4.8)
Blue Sky Bio	One Stage Narrow Platform (4.1), One Stage Regular Platform (4.8)
Clinical House Europe GmbH	PerioType TG (4.8)
Cowellmedi Corporation, Ltd	Atlas Internal Regular, Atlas Internal Taper (All 4.8)
CSM	Apolonia Internal Straight Type, Apolonia Internal Tapered Type (All 4.8)
Dental Ratio Systems (DRS) GmbH	Oktagon (4.8)
Dentis	Internal
Dr. Ihde Dental GmbH	Allfit OSS, Allfit Osseopore, Allfit SSO, Allfit STC, Allfit STI, Allfit STI RT, Allfit STO, Allfit TPG (All 4.8)
Euroteknika	Aesthetica (4.2), Aesthetica (4.8), Aesthetica+ (NP-4.2), Aesthetica+ (RP-4.8)
Hi-Tec Implants	Tapered Tite Fit, Tite Fit (All 4.8)
Implant Direct LLC	SwishPlant (4.8)
Intra-Lock System International	Standard Platform Conic Body (4.3), Standard Platform Straight Body (4.3), Wide Platform Conic Body (5.5), Wide Platform Straight Body (5.5)
Keystone Dental (Lifecore)	Stage-1 RDS
MegaGen	ExFeel Internal Straight (4.8), ExFeel Internal Taper (4.8)
Microtech srl	Microfit (4.8)
MIS Implant Technologies	MISTRAL (4.8)
Neobiotech Co., Ltd	Neoplant Internal (4.8)
Osstem Implant Co., Ltd (Hiossen)	SS I Regular, SS II Regular, SS III Regular (All 4.8)
Renew Biocare AG	Re-Fit TG Octa (4.8), Re-Fit TG ST (4.1)
Southern Implants	ITC, ITSA, ITT, Max IT 7, OCT ITSTA, OCT-XT ITST 12d Co-Axis (All 4.8)
Sterngold Dental, LLC	Stern IC (4.8)
Straumann AG	Hollow Cylinder 15 Degree Angled Esthetic Plus 3.5 (4.8), Hollow Cylinder Esthetic Plus 3.5 (4.8), Standard (RN), Standard Plus (RN), Taper Effect (RN)
Sweden & Martina SPA	Stark-D (4.8)
Sybron Implant Solutions	Endopore One Stage, SybronPRO TL, SybronPRO XRT (All 4.8)
Trinon Titanium	QK
UFIT Implant Inc.	ST1, VT1 (All 4.8)
Zimmer Dental	SwissPlus Straight, SwissPlus Tapered (All 4.8)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.



INTERNAL OCTAGON IMPLANT CHART

6.5mmD INTERNAL OCTAGON

IMPLANT COMPANY	IMPLANTS
AS Technology (Titanium Fix)	
B&W s.r.l.	
Biomain AB	
Biomet 3i	TG OSSEOTITE (6.5)
Blue Sky Bio	One Stage Wide Platform (6.5)
Clinical House Europe GmbH	
Cowellmedi Corporation, Ltd	
CSM	
Dental Ratio Systems (DRS) GmbH	
Dentis	
Dr. Ihde Dental GmbH	Allfit STW (6.5)
Euroteknika	Aesthetica (6.5), Aesthetica+ (WP-6.5)
Hi-Tec Implants	
Implant Direct LLC	SwishPlant (6.5)
Intra-Lock System International	
Keystone Dental (Lifecore)	Stage-1 WDS
MegaGen	
Microtech srl	
MIS Implant Technologies	
Neobiotech Co., Ltd	
Osstem Implant Co., Ltd (Hiossen)	
Renew Biocare AG	
Southern Implants	ITC 6 (6.5), ITT 6 (6.5), Max IT 8 (6.5), Max IT 9 (6.5,7.5), OCT ITST6 (6.5)
Sterngold Dental, LLC	
Straumann AG	Standard (WN), Standard Plus (WN), Taper Effect (WN),
Sweden & Martina SPA	
Sybron Implant Solutions	Sybron PRO TL (6.5), SybronPRO XRT (6.5)
Trinon Titanium	
UFIT Implant Inc.	ST1 (6.5), VT1 (6.5)
Zimmer Dental	

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.



INTERNAL 4SIDE IMPLANT CHART

REGULAR PLATFORM (RP) INTERNAL 4SIDE

IMPLANT COMPANY	IMPLANTS
Blue Sky Bio	Quattro Regular Platform (4.1, 4.8, 5.6)
Straumann	Bone Level Implant Regular Crossfit (4.1, 4.8), Bone Level Tapered Implant Regular Crossfit (4.1, 4.8)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.