



INTERNAL HEXAGON IMPLANT CHART

3.5mmD INTERNAL HEXAGON

IMPLANT COMPANY	IMPLANTS
A.B. Dental Devices	I10, I15, I2 Narrow Platform, I2b Narrow Platform, I5 Narrow Platform (All 3.75)
ACE Surgical	infinity Internal Hex (3.7 & 4.1) (3.5mmD Platform)
ADIN Dental Implants System	Swell, Touareg, Touareg-s
Alfa Gate	Bioactive, DFCI (3.5)
Alpha Bio Tec.	ATID, DFI, ICE (3.5), SFB, SPI
Alpha Dent Implants	Alpha Dent Active, Active Bio, Active Conus, Active Plus, Classic, Classic Conus
American Dental Implant Corporation - ADI	Internal Hex Skinny 2.4, Internal Hex Small Diameter, Internal Hex Standard Diameter, Internal Hex Tapered Standard Diameter (All 3.5)
ARDS Implants	Classic, Dental Implant S, Premium, Smart (All 3.75)
BioHorizon	Internal, Single Stage, Tapered Internal, Tapered Internal Plus (All 3.5)
Blue Sky Bio	Internal Hex (3.5)
CORTEX Dental Implants Industries, Ltd.	Classic, Classix (3.75), Dynamic, Dynamix (3.75)
Den Tack	Quad Internal Hex
DENTIN Implants Technologies LTD	Classic, Prestige, Rapid
Dr. Ihde Dental GmbH	Allfit ATI(3.5)
Easy System Implant	Previum-S (4.1)
Glidewell Laboratories	Inclusive Tapered (3.5)
Hi-Tec Implants	Logic Plus, Smooth & Slim Fit, Tapered Self Thread, X6 (All 3.5)
Implant Direct LLC	Legacy 1, Legacy 2, Legacy 3, Legacy 4 (All 3.5)
Implasa-Höchst GmbH	DoubleGotic, Gotic, Nordic (All 3.5)
IQ Implants	Master, Smart (All 3.5)
Maxillent	iRaise, iSure (All 3.5)
Microtech srl	Minivent (3.5)
MIS Implant Technologies	BioCom Internal Hexagon - Standard Platform, Lance Internal Hexagon Standard Platform, Lance Plus - Standard Platform, M4 Standard Platform, Seven - Standard Platform (All 3.75)
Nova-Implants Ltd	Pacific, Spirit (All 3.5)
Park Dental Research Corp.	Star Vent Internal Hex Screw
Sialo Technology Israel Ltd.	DIVA (Dynamic Implant Valve Approach)
Sigma Implants Ltd	Sigma Active, Sigma Classic
TAV Dental	Silhouette, Sirius (All 3.5)
Tecom Implantology (TITANMED S.r.l)	Screw Straight, Screw Tapered, Starfly 0.2 Two Phases (All 3.5)
Tidal Spiral Implant Systems (TSI)	Tidal Screw Implant System TSS, Tidal Spiral Implant System TSI (All 3.5)
Top Implants LTD.	Dolphin, Orca, T-Shark (All 3.5)
TRI Dental Implants GmbH	TRI-Vent (3.5)
Unidental Direct, LLC	IN-SIX
Zimmer Dental	[Paragon] - Bio-Vent SD (3.5), Core-Vent SD (3.5), Micro-Vent SD (3.25), Micro-Vent Tapered SD (3.7) - [Paragon] AND Screw Vent, Tapered Screw Vent, Trabecular Metal (All 3.5)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.





INTERNAL HEXAGON IMPLANT CHART

4.5mmD INTERNAL HEXAGON

IMPLANT COMPANY	IMPLANTS
A.B. Dental Devices	I1 Wide Platform, I2 Wide Platform, I2b Wide Platform, I5 Wide Platform (All 4.5)
ACE Surgical	infinity Internal Hex (4.7 & 5.1) (4.5 mmD Platform)
ADIN Dental Implants System	
Alfa Gate	
Alpha Bio Tec.	
Alpha Dent Implants	
American Dental Implant Corporation - ADI	Internal Hex Tapered Wide Diameter, Internal Hex Wide Diameter (All 4.5)
ARDS Implants	Dental Implant S (4.5)
BioHorizon	Internal, Single Stage, Tapered Internal, Tapered Internal Plus (All 4.5)
Blue Sky Bio	Internal Hex, Molar (All 4.5)
CORTEX Dental Implants Industries, Ltd.	
Den Tack	
DENTIN Implants Technologies LTD	
Dr. Ihde Dental GmbH	
Easy System Implant	Previum-S (5.0)
Glidewell Laboratories	Inclusive Tapered (4.5)
Hi-Tec Implants	Tapered Self Thread (4.5)
Implant Direct LLC	Legacy 1, Legacy 2, Legacy 3, Legacy 4 (All 4.5)
Implasa-Höchst GmbH	
IQ Implants	
Maxillent	
Microtech srl	Minivent (4.7)
MIS Implant Technologies	BioCom Internal Hexagon Wide Platform, Lance Internal Hexagon Wide Platform, Lance Plus Wide Platform, M4 Wide Platform, One Stage Internal Hexagon, Seven Wide Platform (All 4.5)
Nova-Implants Ltd	Pacific (4.5)
Park Dental Research Corp.	
Sialo Technology Israel Ltd.	
Sigma Implants Ltd	
TAV Dental	
Tecom Implantology (TITANMED S.r.l)	Screw Straight (4.5)
Tidal Spiral Implant Systems (TSI)	Tidal Screw Implant System TSS, Tidal Spiral Implant System TSI (All 4.5)
Top Implants LTD.	Dolphin, Orca, T-Shark (All 4.5)
TRI Dental Implants GmbH	
Unidental Direct, LLC	
Zimmer Dental	[Paragon] - Bio-Vent MD (4.5), Core-Vent MD (4.5), Micro-Vent MD (4.25,4.7), Micro-Vent Tapered MD (4.7) AND Screw-Vent, Tapered Screw-Vent (All 4.5)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.



INTERNAL HEXAGON IMPLANT CHART

5.7mmD INTERNAL HEXAGON

IMPLANT COMPANY	IMPLANTS
A.B. Dental Devices	
ACE Surgical	
ADIN Dental Implants System	
Alfa Gate	
Alpha Bio Tec.	
Alpha Dent Implants	
American Dental Implant Corporation - ADI	Internal Hex Tapered XL (5.7)
ARDS Implants	
BioHorizon	Internal, Single Stage, Tapered Internal (All 5.7)
Blue Sky Bio	
CORTEX Dental Implants Industries, Ltd.	
Den Tack	
DENTIN Implants Technologies LTD	
Dr. Ihde Dental GmbH	
Easy System Implant	
Glidewell Laboratories	
Hi-Tec Implants	
Implant Direct LLC	Legacy 1, Legacy 2, Legacy 3, Legacy 4 (All 5.7)
Implasa-Höchst GmbH	
IQ Implants	
Maxillent	
Microtech srl	Minivent (5.7)
MIS Implant Technologies	
Nova-Implants Ltd	
Park Dental Research Corp.	
Sialo Technology Israel Ltd.	
Sigma Implants Ltd	
TAV Dental	
Tecom Implantology (TITANMED S.r.l)	
Tidal Spiral Implant Systems (TSI)	
Top Implants LTD.	
TRI Dental Implants GmbH	
Unidental Direct, LLC	
Zimmer Dental	[Paragon] - Micro-Vent Tapered WD (5.7) AND AdVent, Tapered Screw-Vent (5.7)

*Documented information has been acquired from 3rd party sources and has not been verified by Buckeye Medical Technologies, LLC.