


Anatotemp[®]

Anatomic Dental Implant
Healing Abutments

"An off the shelf solution for creating
ideal emergence profile."



Easy Online Ordering

Designed for Many Dental
Implant Connections

The Anatotemp anatomic healing abutment aims to replace conventional, round healing abutments. All Anatotemps possess interproximal concavities against which bone or bone grafting material can be placed after the Anatotemp has been completely seated and secured. The smooth medical-grade, plastic healing abutment is held securely by an anti-rotational connection and abutment screw and aids in creating an esthetic emergence through the gingiva during the healing period. Anatotemp aims to be the new standard in dental implant healing abutments.

Six Anatomic Shapes



Numerous Connections



www.anatotemp.com



Buckeye Medical Technologies, LLC.
405 Niles Cortland Rd. SE, Warren, OH 44484
Products shown are not to scale. The Anatotemp logo,
and all other trademarks are trademarks of
Buckeye Medical Technologies, LLC.

©2017 Buckeye Medical Technologies, LLC.

United States Utility Patent #US 8,382,477
Patent also in Canada, Israel, Australia, and the EU

CHOOSE MY ANATOTEMP

Ordering your Anatotemp anatomic temporary dental implant healing abutment is as easy as 1-2-3!

STEP 1

CONFIRM your implant connection.
e.g. Internal Hexagon.

STEP 2

Depending on the tooth being restored, locate the corresponding Anatotemp shape and platform size.

STEP 3

The corresponding box will contain the Anatotemp part number. Place your order online @ www.anatotemp.com.

INTERNAL HEXAGON CONNECTION



	ANATOTEMP SHAPE					
	Mandibular Molar	Maxillary Molar	Maxillary/Mandibular Premolar Mandibular Canine	Maxillary Canine	Maxillary Lateral/Mandibular Incisor	Maxillary Central
Platform Size						
3.5mmD			1035-04	1035-03	1035-02	1035-01
4.5mmD	1045-06	1045-05	1045-04	1045-03		1045-01
5.7mmD	1057-06	1057-05				

Compatible with: Implant Direct LLC - Legacy™ System

INTERNAL TRI-LOBE CONNECTION



	ANATOTEMP SHAPE					
	Mandibular Molar	Maxillary Molar	Maxillary/Mandibular Premolar Mandibular Canine	Maxillary Canine	Maxillary Lateral/Mandibular Incisor	Maxillary Central
Platform Size						
3.5mmD			2035-04	2035-03	2035-02	2035-01
4.3mmD			2043-04	2043-03		2043-01
5.0mmD	2050-06	2050-05				
6.0mmD	2060-06	2060-05				

Compatible with: Implant Direct LLC - ReActive, RePlant, RePlus

INTERNAL CONICAL HEXAGON CONNECTION



	ANATOTEMP SHAPE					
	Mandibular Molar	Maxillary Molar	Maxillary/Mandibular Premolar Mandibular Canine	Maxillary Canine	Maxillary Lateral/Mandibular Incisor	Maxillary Central
Platform Size						
3.0mmD			1530-04	1530-03	1530-02	1530-01
3.4mmD	1534-06	1534-05	1534-04	1534-03		1534-01

Compatible with: Implant Direct LLC - InterActive™ System

* Anatotemp Platform 3.0mm corresponds to NobelActive NP Platform

* Anatotemp Platform 3.4mm corresponds to NobelActive RP Platform

CHOOSE MY ANATOTEMP

Ordering your Anatotemp anatomic temporary dental implant healing abutment is as easy as 1-2-3!

STEP 1

CONFIRM your implant connection.
e.g. Internal Hexagon.

STEP 2

Depending on the tooth being restored, locate the corresponding Anatotemp shape and platform size.

STEP 3

The corresponding box will contain the Anatotemp part number. Place your order online @ www.anatotemp.com.

INTERNAL 4SIDE CONNECTION



		ANATOTEMP SHAPE					
		Mandibular Molar	Maxillary Molar	Maxillary/Mandibular Premolar Mandibular Canine	Maxillary Canine	Maxillary Lateral/Mandibular Incisor	Maxillary Central
Platform Size	Implant Diameter						
RP	4.1mmD			4041-04	4041-03	4041-02	4041-01
RP	4.8mmD	4048-06	4048-05	4041-04	4041-03		4041-01
RP	5.6mmD	4048-06	4048-05				

Compatible with: Blue Sky Bio LLC - Quattro
* Anatotemp Platform RP corresponds to Straumann RP Platform

INTERNAL OCTAGON CONNECTION



		ANATOTEMP SHAPE					
		Mandibular Molar	Maxillary Molar	Maxillary/Mandibular Premolar Mandibular Canine	Maxillary Canine	Maxillary Lateral/Mandibular Incisor	Maxillary Central
Platform Size							
	4.8mmD	3048-06	3048-05	3048-04	3048-03		3048-01
	6.5mmD	3065-06	3065-05				

Compatible with: Implant Direct LLC - SwishPlant